Gebrauchsanweisung	6	DEU
Operating instructions	11	ENG
Mode d'emploi	16	FRA
Istruzioni d'uso	22	ITA
使用説明書	27	JPN

D^{med}® **halux** LED 20-3 P

D^{med}® **halux** LED 20-1 P

Untersuchungsleuchte
Examination light
Lampe d'examens
Apparecchio per visita
検査用ライト

SYMBOLE; SYMBOLS; SYMBOLES ; SIMBOLI; 記号の意味

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The warning symbol indicates all instructions that are important for safety. Failure to comply with them can lead to injury or damage to the luminaire or the equipment. In combination with the following signal words the warning symbol means:

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

CONTENTS

1	SAFETY INSTRUCTIONS	11
1.1	Intended use	11
1.2	User profiles	11
1.3	Safety instructions	11
1.4	Warning levels	11
2	SCOPE OF SUPPLY	11
2.1	halux LED 20-3/-1 P FX	11
2.2	halux LED 20-3/-1 P SX	11
3	MOUNTING	12
4	OPERATION	12
4.1	halux LED 20-3 P FX / SX	12
4.2	halux LED 20-1 P FX / SX	12
5	DISINFECTION CLEANING	12
6	SAFETY INSPECTIONS	13
7	DEMOUNTING	13
7.1	Disposal	13
8	ACCESSORIES	13
9	ADDITIONAL INSTRUCTIONS	13
10	TROUBLESHOOTING	14
11	TECHNICAL DATA	14
12	ELECTROMAGNETIC COMPATIBILITY (EMC)	15

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1 SAFETY INSTRUCTIONS

1.1 Intended use

The D^{med}® halux LED luminaire is an examination light. Its intended use is to illuminate the body of the patient to support diagnosis or treatment. An interruption of the diagnosis or treatment caused by a light failure is possible at any time without risk for the patient. The luminaire is not intended for use in operating rooms.

1.2 User profiles

Medical professional

A medical professional is any person who has completed medical training and is working in his or her professional field.

Cleaning professional

A cleaning professional, as referred to here, must be familiar with national and workplace-related hygiene regulations.

Electrician

An electrician, as referred to here, must have been trained in the fields of electronics and electrical engineering and be familiar with the relevant standards and regulations.

Qualified professional

A qualified professional, as referred to here, must be capable of mounting and demounting the luminaire thanks to professional training, knowledge and experience and knowledge of the relevant regulations.

1.3 Safety instructions

- ▶ Operation by a medical professional
- ▶ The instructions form part of the product. They must be retained and made accessible to all subsequent users.
- ▶ Any work on the luminaire (including repairs) must be carried out by qualified electricians only. Mounting must be performed by a qualified professional only.
- ▶ The luminaire must not be altered or manipulated in any way. Only approved original parts must be used. Any use other than the intended use with original parts may give rise to different technical values and life-threatening danger.
- ▶ Operation in potentially explosive areas is prohibited. The luminaire power supply is a potential ignition source.
- ▶ The luminaire must only be operated in dry, dust-free rooms.
- ▶ The luminaire must not be left switched on without supervision.
- ▶ Do not use any luminaire that is damaged. Defective cables are also potential hazards. Do not place the cable close to any heat source or on sharp edges.
- ▶ Luminaires which are equipped with a hospital grade plug (UL 817, only US- and Canada version), the earth connector is used as a functional earth within the luminaire.
- ▶ Damage to the eyes. Never look directly into the light source.
- ▶ Replace any glass that is damaged before operating the luminaire again.
- ▶ Never place additional loads on the luminaire head or the arm.

- ▶ The luminaire must not be covered by a cloth or any similar item while it is in operation.
- ▶ The ventilation openings (where they exist) must be unobstructed whenever the luminaire is in operation.
- ▶ The luminaire must not be operated near to external heat sources that exceed the maximum ambient temperature for the luminaire.
- ▶ The luminaire must not be used in environmental conditions other than those for which it is intended.
- ▶ The luminaire must only be used for the intended use described in this document.
- ▶ The manufacturer cannot be held responsible for any injury or damage that is caused as a result of any use other than the intended use or of any failure to comply with safety instructions and warnings.

1.4 Warning levels

 DANGER
Indication of hazards that can lead to death or serious injury if instructions are disregarded.
 WARNING
Indication of hazards that can lead to injury if instructions are disregarded.
CAUTION
Indication of hazards that can lead to damage to property if instructions are disregarded.

2 SCOPE OF SUPPLY

2.1 halux LED 20-3/-1 P FX

The following are included in the scope of supply:

- 1x Luminaire
- 1x Power cable

2.2 halux LED 20-3/-1 P SX

The following are included in the scope of supply:

- 1x Luminaire
- 1x Power cable

3 MOUNTING

- ▶ The luminaires are equipped with an adapter pin. This must be used to position the luminaire in one of the accessories described in chapter 8.

Load data:

Bending moment M_B	15 Nm
Vertical weight F_G	15 N

4 OPERATION

⚠ DANGER


Danger of death from electric shock

- ▶ Do not insert any power cable that is damaged.
- ▶ If there is any sign of damage to the power cable, replace it immediately with a new one.
- ▶ The supply voltage and frequency must match the data on the type plate.

⚠ WARNING


Risk of eye damage

- ▶ Never look directly into the light source


- ▶ Insert power cable into the luminaire
- ▶ Connect the power cable to the mains
- ▶ Perform a functional test before every use: all LEDs in the light cone must illuminate.


4.1 halux LED 20-3 P FX / SX

- ▶ Switch the LED module on/off

- ▶ Set the desired illuminance in 20% dimmer steps

- ▶ Set the desired color temperature according to the following table:

COLOR	COL. TEMP.	RECOMMENDATION
Blue (cold white)	4400K	Ideal for identifying individual vessels for wound follow-up treatment and small procedures
White (neutral white)	3800K	Ideal for general examinations
Red (warm white)	3300K	Ideal for various types of skin examinations

4.2 halux LED 20-1 P FX / SX

- ▶ Switch the luminaire on/off

5 DISINFECTION CLEANING

⚠ DANGER

Danger of death from electric shock

- ▶ Switch off the power and secure against unintended switching on before disinfectant cleaning.

CAUTION

Damage to property caused by incorrect cleaning

- ▶ Use only disinfectants that will not affect the functional capability of the light.
- ▶ No cleaning agents containing solvents, chlorine or abrasive products must be used. Those agents can cause damage to the plastic parts.
- ▶ The cleaning and disinfection agents must be suitable and approved for plastics like PC, PMMA, PA and ABS.
- ▶ Concentrated disinfectants may damage the light.
- ▶ Pay attention to the specifications from the data sheet of the agent for concentration and effectiveness time.
- ▶ Unsuitable cloths can cause scratches.

CAUTION

Dirt decreases luminosity

- ▶ Clean regularly to keep the cover clean.
- ▶ Only disinfection by wiping is permitted.

- ▶ Clean the PC cover with a suitable cleaning cloth and a suitable cleaning agent.

CAUTION

To minimize the risk of disease transmission, in addition to complying with this user manual, you must also comply with the applicable occupational health and safety regulations and the requirements of national bodies with responsibilities for hygiene and disinfection.

6 SAFETY INSPECTIONS

DANGER

Danger of death from electric shock

- ▶ Remove the plug from the mains and put the switch in the "OFF" position.
- ▶ The power supply cable must be checked for damage at least once per year.

CAUTION

- ▶ Maintenance and repairs must only be carried out by qualified electricians.
- ▶ The corresponding user profile is described in Chapter 1, Safety instructions.

7 DEMOUNTING

DANGER

Danger of death due to electric shock

- ▶ Disconnect from the power supply and secure against unintended switching on before disassembling.

7.1 Disposal

Do not dispose of the luminaire in household refuse. Dispose of the luminaire at a disposal point in accordance with local regulations or take them to a dealer that provides an appropriate disposal service. Cut off the cable at the housing.

The products listed above are more than 95% recyclable. The luminaires have been constructed to be compatible with recycling so that a high proportion of the materials used in these products can be recycled or converted into energy after the end of life cycle. They contain no materials that are dangerous or that need to be monitored.

8 ACCESSORIES

Table clamp (Art. no.: D14.228.000 - black)

Universal mount (Art. no.: D13.430.000 - pure white)


Rail clamp (Art. no.: D13.269.000)

Wall bracket (Art. no.: D13.231.000 - pure white)

Roller stand RL (Art. no.: D14.022.000)

9 ADDITIONAL INSTRUCTIONS

The luminaire itself is maintenance-free.

You can obtain additional documents about this product from the manufacturer upon request. These luminaires do not affect any other devices. To save energy, the luminaire should only be switched on when actually needed.

10 TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profiles
The luminaire does not come on	Contact problem	Switch on again	All
The luminaire does not come on	No mains voltage	Check mains voltage, check all connections	Electrician
The luminaire does not come on	LED module defective	Contact manufacturer support	By manufacturer service only

11 TECHNICAL DATA

Electrical data:

Rated input voltage	100-240V or 120V or 220-240V see luminaire type plate
Frequency range	50/60Hz
Power consumption (20-3) (20-1)	15-22VA (12-18W) 9VA / 15VA (8.5W / 9.5W)
Input current (20-3) / (20-1)	0.09A-0.15A / 0.07A

Photometric values:

Central illuminance E_c at 0.5m (1.64 feet) distance	40'000 lx *
Light field diameter d_{10} at 0.5m (1.64 feet) distance	$\varnothing = 18 \text{ cm (7.1")} *$
Light field diameter d_{50} at 0.5m (1.64 feet) distance	$\varnothing = 9 \text{ cm (3.54")} *$
Color temperature (20-3) (20-1)	4400K / 3800K / 3300K * 4400K *
Color rendering index R_a	>93
Color rendering index R_9	>90
Total irradiance E_e at maximum intensity	<180 W/m ² * -10% / +20% tolerance

Ambient conditions for transport, storage and operation:

Ambient temperature (storage and transport)	-20°C to +70°C (-4°F to +158°F)
Ambient temperature (operation)	+10°C to +35°C (+50°F to +95°F)
Relative humidity (non-condensing) (storage and transport)	max. 90%
Relative humidity (non-condensing) (operation)	max. 75%

Weight:

halux LED 20-3/-1 P FX	1.4kg (3.09 lb)
halux LED 20-3/-1 P SX	1.4kg (3.09 lb)

Operating mode:

Operating mode	Continuous operation
----------------	----------------------

Classification:

halux LED 20-1 P FX / SX	Protection class II
Degree of protection as per IEC 60529	IP 20
Classification according to Directive 93/42/EEC - Annex IX (medical device classes)	Class I
Electrical safety test and EMC according to:	EN/IEC 60601-1 EN/IEC 60601-2-41 EN/IEC 60601-1-2
GMDN Code	12276

Service life of light source:

LED life cycle	30'000h (L70/B50)
----------------	-------------------

12 ELECTROMAGNETIC COMPATIBILITY (EMC)

Electromagnetic emission guidelines		
The medical device is intended for operation in an electromagnetic environment such as the one specified below. The user must ensure that it is operated in such an environment.		
Emissions	Compliance	Electromagnetic environment
High-frequency emissions (CISPR 11)	Group 1	The medical device uses HF energy exclusively for its internal function. This means that its HF emissions are very low, and it is unlikely that adjacent electronic devices will experience interference.
High-frequency emissions (CISPR 11)	Class B	The medical device is intended for use in all facilities, including residential buildings and facilities that are directly connected (without a transformer) to the same low voltage network as residential buildings.
Emissions from harmonics (IEC 61000-3-2)	Class C	
Emissions from voltage fluctuations/flicker (IEC 61000-3-3)	Compliant	

Electromagnetic interference immunity guidelines			
The medical device is intended for operation in an electromagnetic environment such as the one specified below. The user must ensure that it is operated in such an environment.			
Immunity to	IEC 60601-1-2 testing level	Compliance level of the medical device	Electromagnetic environment
Electrostatic discharge (ESD) (IEC 61000-4-2)	Contact discharge: ± 6 kV Air discharge: ± 8 kV	± 6 kV ± 15 kV	Floors made of timber, concrete or ceramics tiles are preferred. Where synthetic floor covering is used the relative humidity should be at least 30%.
Electrical fast transients/bursts (IEC 61000-4-4)	Power supply cables: ± 2 kV Longer input and output cables: ± 1 kV	± 2 kV Not applicable	The quality of the supply voltage should match that of a typical business or hospital environment.
Surge voltages/surges (IEC 61000-4-5)	± 1 kV phase-to-phase voltage ± 2 kV phase-to-ground voltage	± 1 kV n/a	The quality of the supply voltage should match that of a typical business or hospital environment.
Magnetic field in the power supply frequency (50/60 Hz) (IEC 61000-4-8)	3 A/m	100 A/m	In proximity of this medical device, do not operate devices with unusually strong network-frequency magnetic fields (transformer stations etc.).
Voltage sags and short interruptions to the supply voltage (IEC 61000-4-11)	Sag >95 %, 0.5 cycles Sag 60 %, 5 cycles Sag 30 %, 25 cycles Sag >95 %, 5 seconds	Sag >95 %, 0.5 cycles Sag 60 %, 5 cycles Sag 30 %, 25 cycles Sag >95 %, 5 seconds	The quality of the supply voltage should match that of a typical business or hospital environment. If the user requires continued function during any interruption of the energy supply system, we recommend that the medical device be powered by an uninterrupted power supply or a battery.
Radiated HF disturbances (IEC 61000-4-3)	3 V/m 80 MHz - 2.5 GHz	10 V/m	Recommended minimum distance of portable and mobile radio devices of radiated power PEIRP to the medical device, including its cables: $d = 0.35\sqrt{P}$
Conducted HF disturbances (IEC 61000-4-6)	3 $V_{\text{rms-value}}$ 150 kHz - 80 MHz:	10 $V_{\text{rms-value}}$	Recommended minimum distance of portable and mobile radio devices of radiated power PEIRP to the medical device, including its cables: 80 MHz - 800 MHz: $d = 0.35\sqrt{P}$ 800 MHz - 2.5 GHz: $d = 0.7\sqrt{P}$
d = recommended safe distance [m], P = Rated Power of the sender [W]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level for all frequencies.			
Interference may occur in the vicinity of equipment marked with: 			

Recommended safe distances from portable and mobile HF communications equipment		
Rated Power of the transmitter [W]	150 kHz - 800 MHz	800 MHz - 2.5 GHz
0.01	0.035 m (0.11 ft)	0.07 m (0.23 ft)
0.1	0.11 m (0.36 ft)	0.22 m (0.72 ft)
1	0.35 m (1.15 ft)	0.7 m (2.30 ft)
10	1.11 m (3.64 ft)	2.21 m (7.25 ft)
50	2.47 m (8.10 ft)	4.95 m (16.24 ft)
100	3.5 m (11.48 ft)	7 m (22.97 ft)

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Derungs Licht AG
Herbert Waldmann GmbH & Co. KG